Austausch-/Reparaturanleitungen * Exchange and repair procedures**

Wechsel des Membranbeutels (http://ifaerg.de/server/downloads/bladder_exchange.mp4):

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Bild 4 (a, b):



Bild 5 (a, b):



Exchange of membrane bladder (http://ifaerg.de/server/downloads/bladder_exchange.mp4):

- Loosen the plug-in flange connection between the bladder ring and the adapter unit (red) by turning and careful levering off (screwdriver) against the angled hose connection pieces. See Fig. 1: a, b
- Strong pushing out of the bladder ring from the adapter unit as well as from the distribution system with sensor (folding of the silicone tubing for the passage through the ring). See Fig. 2: a, b
- **Put O-Ring with sufficient sliding lubricant into the corresponding nut inside of the bag ring**
- Insertion of the distribution system with sensor into the new bladder. See Fig. 3: a, b
Hint: blow some air into the bladder
Check: O-ring inserted in the bag ring?
Check: rubberized membrane side (no vapour transmission) towards the torso central plate
- Put sliding lubricant inside the bag ring peripherally over the O-ring
- Put the bladder ring on the adapter unit by turning and lightly pushing into correct position. See Fig. 4: a, b
- **Note**: only grasp the rubber ring, **do not press on the edge of the ring (membrane damage!)**; see Fig. 5a
- Check correct position of distribution system with sensor inside the bladder; see Fig. 5b
- When inserting the cover shells, **avoid squeezing of the bladder in the groove/nut!**
- Filling and commissioning

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
Exchange of HPK-unit:

- The new HPK can be directly replaced just by setting into the corresponding frame in the torso.
 - Remove bladder-support-plate (dark grey, 6 screws) and loosen the red bladder-adapter-unit (2 screws, bladder ring remains)
 - Pull off both HPK-tubes, disconnect plug and remove HPK
 - Insert new HPK correspondingly (tubes connections must not be changed!), stick necessary distance rubbers, connect plug and place cable (see picture)
 - Fix Bladder support plate and fix bladder adapter unit
- Check finally the correct filling by listening to the sound of the pump (disconnect the other unit because both pumps are controlled simultaneously). During operation the sound should be constant (monotonic), not amphoric or vibrating (when unfilled) and not slurping (water level is too low).
- If there is a damage of the SCU connection cable (check with SPY-adapter) please replace as follows before replacing the HPK: cut off the 5-pin-SCU3 plug and loosen the grommet in the plexi-plate of the torso. Take off the broken cable and replace it with the new. If the grommet has also broken, replace it with the spare part. Adjust the cable length inside (to reach HPK-plug) and outside the torso and fix/tighten the grommet. Finally solder the 5-pin-plug (spare part) to the free end of the cable as follows: brown (pump Vcc) to pin #4, green to pin #5 (pump gnd), yellow to pin #2 (heater Vcc), white to pin #1 (heater gnd), pin #3 is unused.
- When operating the Torso check periodically for sufficient water level. Use only the filling nozzle close to the center plate (see picture, that with the pit) and connect the others to ensure a complete air removing. Fill in until water level at all three nozzles is the same.

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Filling and commissioning of the sweating units:

- At regular intervals, i. daily at 8h continuous operation and twice a week at about 4h operation per day and at uneven pumping noise, check the water level for completeness and refill if necessary.
- For this, use the tube connector that is located closest to the torso central plate (see arrow in figure). Also open the other sockets by attaching the fittings so that air can escape
- Add water until all 3 connections show the same water level.
- Finally check filling by short test operation. In doing so, the pump (possibly disconnecting another unit on the SCU, as both pumps are controlled in parallel) generates a monotonous, even running noise that should not sound amphoric or vibrating or slurping. Otherwise refill again.
- **Pump rate and corresponding filling level are monitored on the process screen of SCU: MIN-low rate (by air!), ON-correct rate, MAX-blocking**
- The problem of low water level does not occur when the associated unit is at the top (usually the ventral side) because HPK is the lowest point in the system. However, there must still be enough water in the subsystem so that the return flow can take place via the red adapter unit. In order to check complete filling, it is advisable to turn upwards to the associated subunit and to ensure a uniform pump sound. If the pump noise remains steady when turning back, filling is sufficient.